



Pulsed Light Decontamination PASS BOX

Designed to decontaminate material before entering a clean zone

FEATURES

- Inlet or outlet Pass Box
- Interlocked door
- Discontinuous system
- Self-decontamination of the chamber
- · Air-tight, overpressure on the cleaner side
- No shadow area when placing the object to be treated
- Several sizes of machine available, adapted to the objects to be treated

BENEFITS

- Efficiency: 6 log on Bacillus pumilus
- Reproducible decontamination process
- Treatment < 1 min
- Low energy consumption
- Free from chemical and residue, no radiation

Discover the Claranor Pharma Pulsed light range

